

CARING FOR YOUR EMPTIER

(mini emptier Lite, R-Net)

Your new emptier is very reliable and requires little maintenance. Clean water is what we recommend for a daily flushing. You may use a germicide or vinegar as long as you flush with water last. Vinegar will damage the plunger pad if left in the emptier. Chlorine bleach should never be used.

The only equipment needed is a squeeze bottle with the same size hose as the leg-bag. The following should be incorporated into your end of day routine.

- 1) Empty leg-bag, turn off or clamp manual valve.
- 2) With emptier still on, have your attendant pull the leg-bag hose off the emptier while at the same time holding down the emptier with other hand.
- 3) Replace the leg-bag hose (on the emptier) with the hose on the cleaner bottle, invert the bottle and flush the valve with clean water, for a few of seconds.

On the rare occasion that your emptier may drip or leak, it will require some minor disassembling and cleaning. Usually this is caused by calcium build-up or pubic hairs interfering with the emptier seal. The following steps will need to be taken.

- 1) Flush emptier as outlined above.
- 2) Remove from mounting location or from hose if mounted inline.
- 3) Follow the steps as shown in the attached Mini Lite Cleaning Sheets.

TROUBLE SHOOTING

1. -Check that all connectors are plugged in properly.
-Check that all wires are intact and not broken or bare.
2. -Check that the R-Net plug is firmly seated.
3. -Check the fuse. (1A 250V fast blow). Located inside the R-Net adapter.
4. -If the fuse blows for an unknown reason have the system checked by a qualified electrical technician before operating again.
5. -If the system drains unexpectedly, or the emptier will not shut off, the switch(s) harness has likely been pinched and is shorting. Disconnect the R-Net plug and have the system checked by a qualified electrical technician before operating again.
7. -Contact:
Assistive Designs
P.O. Box 20015
55 Muskoka Rd. 118 West
Bracebridge, Ontario
Canada, P1L 2H6
fax: (705) 646-8845
e-mail: emptiers@assistivedesigns.com
www.assistivedesigns.com

WARRANTY - Assistive Designs warrants each of its products to be free from defects in material and/or workmanship for a period of twelve (12) months from the date of purchase by the customer. This warranty is expressly in lieu of all other warranties, expressed or implied. This warranty shall not apply to Assistive Designs equipment which shall have been repaired or altered by unauthorized service personnel, or which has been subjected to misuse, damage, or operated in a way other than its intended use. This warranty does not cover damage in transit. ASSISTIVE DESIGNS, ITS EMPLOYEES AND OFFICERS SHALL NOT BE LIABLE FOR ANY SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES, OR RESULTANT COSTS HOWEVER CAUSED.